



DEKO

DC2200H | Drying cabinet

FRANKE

Anaesthetic hoses

- 6+6 pcs cassettes and
- holder shelf for 36 long AN-hoses

Five programmable cycles

Time and temperature adjustment

Main switch

Hepa-filter


The sound and heat insulated double walls of the cabinet are made of AISI304 stainless steel

DC2200H is available as a single door or pass-through model, left or right-hinged


Standard equipment:
8 rails and
3 shelves



DC2200H Drying cabinet is designed especially for drying surgical instruments, anaesthetic hoses & equipment, as well as other re-usable items used in patient health care.



DC2200H Drying cabinet



Electricity: 400V 3N, 16 A, 5,6 kW

Air exhaust: connection 100 mm; 100 m³/h

Wall opening on double door version: W850 x H2020 mm

Weight: 185 kg





DEKO

Franke Medical Oy
Vartiokuja 1
FI-76850 Naarajärvi
Finland
Tel. +358 15 34 111
www.frankemedical.com
Email: ws-medical.fi@franke.com

FRANKE